

TEST PROCEDURE STEPS

<p>① Wash</p>	<p>② Peel</p>	<p>③ Place</p>	<p>④ Insert tip</p> <p>2 to 4cm</p>
<p>Wash your hand before start test.</p>	<p>Peel off the aluminum foil on the extraction tube.</p>	<p>Place the tube in the work station.</p>	<p>Gently insert the swab into one nostril. Insert the swab tip 2-4 cm (1-2 cm in children).</p>
<p>⑤ 10 sec. x5</p> <p>1st nostril</p>	<p>⑥ 10 sec. x5</p> <p>2nd nostril</p>	<p>⑦ Rotate 10 sec.</p>	<p>⑧ Remove</p> <p>Squeeze → Squeeze</p>
<p>Rotate the swab along the nasal mucosa 5 times within 7-10 seconds.</p>	<p>Repeat the process with the same swab in the other nostril.</p>	<p>Soak the sampling swab into the extraction reagent. Rotate the swab and press for about 10 seconds.</p>	<p>Remove the swab while squeezing the tube side to expel any liquid.</p>
<p>⑨</p>	<p>⑩ Remove</p>	<p>⑪ 2 Drops</p>	<p>⑫</p> <p>COVID-19 Antigen</p>
<p>Place the extraction tube cover firmly on the extraction tube.</p>	<p>Open the aluminum foil pouch of the test cassette.</p>	<p>Add 2 drops of the treated sample into the sample well.</p>	<p>Read the test result after adding the sample for 10 minutes.</p>

⑬ Read the results

Please share your test result with your healthcare provider and carefully follow your local COVID guidelines/requirements

POSITIVE RESULT

One coloured line should be in the control region (C) and another coloured line should be in the Test region (T).

Note: The intensity of the colour in the test line region (T) will vary based on the amount of SARS-CoV-2 antigen present in the sample. So any shade of colour in the test region (T) should be considered positive.

A positive result means it is very likely you have COVID-19, but the positive samples should be confirmed to reflect this. Immediately go into self-isolation in accordance with the local guidelines and immediately contact your general practitioner/doctor or the local health department in accordance with the instructions of your local authorities. Your test result will be checked by a PCR confirmation test and you will be explained the next steps.



NEGATIVE RESULT

No apparent coloured line appears in the test line region (T). You are unlikely to have COVID-19. However, it is possible for this test to give a negative result that is incorrect (a false negative) in some people with COVID-19. This means you could possibly still have COVID-19 even though the test is negative.

If you experience symptoms such as headaches, migraines, fever, loss of sense of smell or taste, contact the nearest medical facility according to the rules of your local authority. In addition, you can repeat the test with a new test kit. In case of suspicion repeat the test after 1-2 days, as the corona virus cannot be precisely detected in all phases of an infection. Even with a negative test result, distance yourself and hygiene rules must be observed, migration/traveling, attending events, etc. you should follow your local COVID guidelines/requirements.




INVALID RESULT (Test did not work)

Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test or contact our COVID-19 test centre.



⑭



Put all used test kit components into the biohazard waste bag and throw away it in the trash.