en

## MONOVISC®

### Lightly Cross-linked High Molecular Weight Sodium Hyaluronate

This instructions for use is intended exclusively for distribution outside of the USA.

In an effort to best meet the needs of our customers and minimize waste, these Instructions for Use are being provided in electronic format. This document is subject to change; the most current version of this IFU is available online. If unsure if using the latest revision, please reprint the IFU at [www.anikaifu.com](http://www.anikaifu.com). If a paper copy is preferred it may be requested, free of charge, by contacting Anika Therapeutics, Inc. at [www.anikaifu.com](http://www.anikaifu.com). The onus resides with the user to ensure that the most up-to-date IFU is used.

Products may not be available in all markets. Product availability is subject to the regulatory clearance in individual markets. Please contact your local representative if you have questions about product availability in your area.

#### **DESCRIPTION:**

MONOVISC® is a sterile, non-pyrogenic, sodium hyaluronate solution, lightly cross-linked with a proprietary chemical cross-linker. MONOVISC contains 22 mg/mL of lightly cross-linked sodium hyaluronate (NaHA) dissolved in phosphate buffered saline at physiological osmolality. MONOVISC is manufactured from ultra pure, high molecular weight sodium hyaluronate produced by bacterial fermentation. Hyaluronic acid is a natural complex polysaccharide of the glycosaminoglycan family.

#### **CHARACTERISTICS:**

Sodium Hyaluronate is a high molecular weight polysaccharide composed of sodium glucuronate and N- acetylglucosamine. The sodium hyaluronate in MONOVISC is derived from bacterial fermentation. Hyaluronic acid is ubiquitously distributed throughout the tissue of the body and is present in high concentrations in such tissues as vitreous humor, synovial fluid, umbilical cord and dermis. Sodium hyaluronate functions as a tissue lubricant and is thought to play an important role in modulating the interactions between adjacent tissues. It can also act as a viscoelastic support maintaining a separation between tissues. Different sodium hyaluronate preparations may have different molecular weights, but have the same chemical structure. The MONOVISC cross-linked sodium hyaluronate injection is biocompatible, non-inflammatory and non-pyrogenic. Sodium hyaluronate preparations have been shown to be well tolerated in osteoarthritic synovial joints such as the knee.

#### **APPLICATION:**

MONOVISC Cross-Linked Sodium Hyaluronate Injection is a single, intra-articular injection of cross-linked sodium hyaluronate designed to treat the symptoms of osteoarthritis.

#### **INDICATIONS:**

MONOVISC is indicated as a viscoelastic supplement or a replacement for synovial fluid in human knee joints. MONOVISC is well suited for treatment of the symptoms of mild to moderate osteoarthritis of the human knee joints for patients who have failed to respond adequately to conservative non-pharmacologic therapy and simple analgesics. In clinical studies, Monovisc has been proven to have a duration of effect of at least 6 months and a second injection 6 months after the first injection was shown to be safe.

#### **DIRECTIONS FOR USE:**

The required amount of MONOVISC is injected through a sterile, disposable, hypodermic needle of suitable gauge into the selected joint space. The sterile needle should be attached to the MONOVISC syringe by a health care professional using a health care facility-approved aseptic technique. Common needle gauges for injections are 18-21 gauge. The final needle selection for any procedure is determined by the physician. The health care provider should ensure proper penetration into joint synovial space prior to injecting MONOVISC.

## **CONTRAINDICATIONS:**

MONOVISC is composed of cross-linked sodium hyaluronate and may contain trace amounts of gram positive bacterial proteins. The following pre-existing conditions may constitute relative or absolute contraindications to the use of MONOVISC:

- Known sensitivity to any of the materials contained in MONOVISC
- Pre-existing infections of the skin region of the intended injection site
- Known infection of the index joint
- Known systemic bleeding disorders
- Weak populations, such as children, pregnant or lactating women

## **PRECAUTIONS:**

- Those precautions normally considered during injection of substances into joints are recommended.
- Only medical professionals trained in accepted injection techniques for delivering agents to joint spaces should inject sodium hyaluronate for this application.
- The amount of MONOVISC® necessary to be injected depends on patient anatomy and needs to be defined by the medical professional performing the procedure. An excess quantity of sodium hyaluronate should not be used and the patient should be monitored closely.
- The synovial space should not be overfilled.
- If pain increases during the injection procedure, the injection should be stopped and the needle withdrawn.
- Transient pain or swelling may occur after the intra-articular (IA) injection.
- As with any invasive joint procedure, it is recommended that patients avoid strenuous or prolonged (i.e., more than one hour) weight-bearing activities such as running or tennis within 48 hours following the intra-articular injection.

## **SSCP**

The summary of the safety and clinical performance (SSCP) is available in EUDAMED (<https://ec.europa.eu/tools/eudamed>) which is linked to the product Basic UDI- DI, 00817337MONOVGU.

## **ADVERSE REACTIONS:**

Adverse reactions to intra-articular injection of Monovisc based on clinical study data tend to be transient in nature and include the following:

- Joint Pain/Arthralgia
- Injection Site Pain
- Swelling
- Stiffness
- Pain/Pain in Extremity
- Effusion
- Injection Site Reaction
- Bruising
- Rash
- Synovitis
- Subcutaneous Nodule
- Cyst
- Redness

## **HOW SUPPLIED:**

MONOVISC is a sterile viscoelastic preparation supplied in a disposable glass syringe delivering 4.0 mL. Each mL of MONOVISC contains 22 mg/mL of lightly cross-linked sodium hyaluronate (NaHA) dissolved in phosphate buffered saline. Note: The contents of the syringe are sterile; however, the product tray is non-sterile.

**FOR INTRA-ARTICULAR USE. STORE AT 2°C to 25°C. PROTECT FROM FREEZING.**

**CAUTION: This device is restricted to sale and use by or under the supervision of a physician.**

This product is for single patient use only and must not be re-sterilized. Reuse of needles or syringes used to inject this product can result in transmission of infectious agents as well as blood-borne pathogens (including HIV and hepatitis), potentially endangering patients and physicians and staff. Used needles or syringes should be discarded after each injection session and not saved for subsequent sessions on the same patient.

**DO NOT USE IF PRODUCT INNER PACKAGING IS OPEN OR DAMAGED.**

To report adverse events or serious incidents or opened/damaged poorly performing products, contact Anika Therapeutics, Inc. at +1-781-457-9000 or [GlobalComplaints@anika.com](mailto:GlobalComplaints@anika.com). In case of any serious incident that has occurred in the EU, Anika Therapeutics should also contact the competent authority of your member state.

	European conformity and notified body number		Batch, LOT number
	Sterilized using steam or dry heat		Catalog Number
	Single barrier system with protective packaging outside		Medical Device
	Single barrier		Do not use if package is damage or opened
	Do not re-use		Importer
	Do not re-sterilize		Prescription only
	Use by		Consult Instructions for use
	Temperature Limit	<a href="http://www.anikaifu.com">www.anikaifu.com</a>	
	Date of manufacture		EU Authorized Representative Anika Therapeutics, S.r.L
	Legal Manufacturer Anika Therapeutics, Inc. 32 Wiggins Avenue Bedford, MA 01730 USA		Corso Stati Uniti, 4/U 35127 Padova (PD) - Italy

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