NOVAMA® MultiCheck Pro+

TOTAL CHOLESTEROL TEST STRIPS

For use with NOVAMA MultiCheck Meter series.

IMPORTANT: Please read this instruction and user's manual before using NOVAMA MultiCheck Pro+ Total Cholesterol Test Strips (Model: BK-C2).

INTENDED USE

- Measuring total cholesterol level in capillary whole blood from fingertip.
- The meter can be used by laypersons or healthcare professionals.
- In vitro diagnostic use only.

PRINCIPLE OF THE TEST

The blood sample is drawn into the strip reaction zone by capillary force. When the filtered blood reacts with the reagent on the electrodes, the meter will detect a current and provide a potential. Then the current will be converted to plasma equivalent total cholesterol value.

CHEMICAL COMPOSITION OF THE STRIP

- Total Cholesterol esterase: ≤ 2.5 Unit
- Total Cholesterol oxidase: ≤ 2.5 Unit
- Peroxidase: ≤ 1 Unit
- Other ingredients (buffer, mediator, surfactant, stabilizer, etc.): ≤ 0.02 mg

STORAGE AND HANDLING

- Store the strips within 4-30°C (39-86°F), use within 3 months after first open.
- Handle the strips with dry and clean hands.
- Recap strip vial immediately after use if applicable.
- Use the strip immediately after taking out from the vial or the foil pack.
- Do not use expired strips. It may cause false result.
- Do not expose the strips to sunlight, heat or humid environment.
- Do not repack, change vial or the foil pack, refrigerate or freeze the strips.
- Do not bend, cut or fold the strips.

AUTO STRIP RECOGNITION

Please refer to User's Manual for more information.

CODING PROCEDURE

- Insert the code strip into the NOVAMA MultiCheck meter.
- Code number and test mode displays on screen.
- Make sure the code on screen, code strip, strip vial label or the foil pack is the same.
- Make sure the test mode on top of the display shows "CHOL".

HOW TO PERFORM THE TEST

- Take a test strip (recap strip vial immediately if applicable) and insert the strip into the meter.
- Check if the code is correct.
- Use the lancing device to prick your finger to create a capillary blood sample.
- Touch the blood sample with strip and make sure blood fills the inlet sufficiently.
- Result will show on the monitor after count down.
- Discard the used strip and lancet to waste container. Please refer to User's Manual for detailed procedure.

IMPORTANT: Please double check the test mode before testing.

TEST RESULTS

- Measuring Range: 100-400 mg/dL (2.59-10.35 mmol/L)
- The Normal Range for Healthy Adult: ≤ 200 mg/dL (5.17 mmol/L). (The range is only a reference.)
- If "Lo" is displayed on the meter, total cholesterol level may be below 100 mg/dL (2.59 mmol/L).
- If "Hi" is displayed on the meter, total cholesterol level may be over 400 mg/dL (10.35 mmol/L).

If the test result is unusually low, high or inconsistent with how you feel, please check the following and repeat the test with a new strip making sure you have followed the instructions carefully:

- Check if strips are not expired.
- Make sure the meter has been coded properly, and the code number is consistent with the strip.
- See if the meter and strips are working properly.
- Make sure that blood sample has filled up the strip inlet sufficiently.

If you still have doubt with the test result, please consult healthcare professional or contact authorized distributor.

IMPORTANT:

- △ Please consult healthcare professional before making medical decision based on the tested results.
- △ Strips from other brands are not compatible with NOVAMA MultiCheck series meters.
- △ POTENTIAL BIOHAZARD. The used test materials could be sources of infection, please beware especially when the system is used by multiple users.

CONTROL SOLUTION TEST

Control solution is used to check the performance of the system. If the control solution test result falls within the range listed on the strip vial label or the foil pack, the system is working properly. Control solution is not included in the system. Please contact an authorized distributor for purchasing.

A CONTROL TEST SHOULD BE PERFORMED:

- When the system does not work properly.
- When the test result is unusual or inconsistent.

OUT OF RANGE CONTROL SOLUTION TEST RESULT MAY BE CAUSED BY:

- Performing the test inaccurately.
- Contaminated or expired control solution.
- System deterioration or malfunction.
- The codes between meter and strip are inconsistent.

IMPORTANT:

- △ Each strip vial or the foil pack might have different control range. Please check every time when you perform a control solution test. Do not use the system if the control solution tests keep producing wrong results. Please contact authorized distributor for assistance.

LIMITATIONS

For getting an accurate test result, please read the following limitation carefully.

- Single use only. Do not reuse the strips.
- Do not perform the test under strong electric field with signal interference.
- Do not perform the test on newborns.
- The test strip can be used at altitude up to 10,000 feet (3048 meters).
- Hematocrit (HCT) should be between 30-50%.
- Use only fresh capillary whole blood from fingertip.
- Do not use serum, plasma or venous whole blood.

ADDITIONAL INFORMATION FOR HEALTHCARE PROFESSIONALS

- Interferences: Genticic acid at therapeutic levels ≥ 2.5 mg/dL significantly affect results.
- In situation of decreased peripheral blood flow false results may occur. Ex: severe dehydration; in shock, or in a hyperosmolar state (with or without ketosis), hypotension.
- Triglyceride level over 1500 mg/dL (16.94 mmol/L) is not recommended for testing with cholesterol test strips.
- Glucose levels up to 500 mg/dL (27.9 mmol/L) do not have significant interference on cholesterol test.
- The test strips are not recommended for critically ill patients.

ACCURACY:

The accuracy of the total cholesterol was performed with capillary finger blood and compared with reference instrument calibrators, which were referenced to the NIST SRM 913 Standard.

At total cholesterol concentration < 150 mg/dL (< 3.89 mmol/L):

Within ± 10 mg/dL (Within ± 0.26 mmol/L)	Within ± 15 mg/dL (Within ± 0.39 mmol/L)	Within ± 20 mg/dL (Within ± 0.52 mmol/L)	Within ± 30 mg/dL (Within ± 0.78 mmol/L)
19/24 (79.2%)	21/24 (87.5%)	23/24 (95.8%)	24/24 (100%)

At total cholesterol concentration ≥ 150 mg/dL (≥ 3.89 mmol/L):

Within $\pm 5\%$	Within $\pm 10\%$	Within $\pm 15\%$	Within 20%
33/99 (33.3%)	58/99 (58.6%)	77/99 (77.8%)	94/99 (94.9%)

PRECISION:

Repeatability:

Mean (mg/dL)	143	185	247
Mean (mmol/L)	3.70	4.80	6.38
SD (mg/dL)	7.0	9.8	10.9
SD (mmol/L)	0.18	0.25	0.28
CV (%)	4.9	5.3	4.4

Intermediate precision:

Mean (mg/dL)	155	247
Mean (mmol/L)	4.01	6.39
CV (%)	4.3	4.7

SYMBOLS

Please refer to User's Manual for more information.

REFERENCE

- Executive Summary of The Third Report of The National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, And Treatment of High Blood Cholesterol In Adults (Adult Treatment Panel III). JAMA. 2001;285:2486-2497.

General Life Biotechnology Co., LTD.
5F., No. 240, Shinshu Rd., Shin Juang Dist., New Taipei City 242, Taiwan

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General Life Biotechnology Co.,LTD.
MDSS GmbH